



Recombinant *Saccharomyces cerevisiae*

Heterogeneous nuclear rnp K-like protein 2 (HEK2)

Product Code	CSB-BP511359SVN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C8Z3W4
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSQFFEAATP VAIPTNNTNG GSSDAGSAAT GGAPVVGTTA QPTINHRLLL SLKEAAKIIG TKGSTISRIR AANSVKIGIS EKVPGCSDRI LSCAGNVINV ANAIGDIVDV LNKRNPNED AAEGEAEHHY YFHFLNHILP APSKDEIRDL QQLEDIGYVR LIVANSHISS IIGKAGATIK SLINKHGVKI VASKDFLPAS DERIIEIQGF PGSITNVLIE ISEILSDVD VRFSTERSYF PHLKKSSGEP TSPSTSSNTR IELKIPELYV GAIIGRGMNR IKNLKTFTKT NIVVERKDDD DKDENFRKFI ITSKFPKNVK LAESMLLNKLN NTEIEKRENY KRKLEAAEGD ATVVTERTSDS ASFLEEKEEP QKNHNDNKEEQ S
Source	Baculovirus
Target Names	HEK2
Protein Names	Recommended name: Heterogeneous nuclear rnp K-like protein 2 Alternative name(s): KH domain-containing protein 1
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.