



Recombinant *Citrobacter rodentium* Transcriptional regulator NanR (nanR)

Product Code	CSB-MP511441DTP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	D2TP58
Product Type	Recombinant Protein
Immunogen Species	<i>Citrobacter rodentium</i> (strain ICC168) (<i>Citrobacter freundii</i> biotype 4280)
Purity	≥85% (SDS-PAGE)
Sequence	MKTFDSQAEN SPAVIGRSLR NRPLARKKLS EMVEEELEQM IRRKEFGEGE QLPSERELMA FFNVGRPSVR EALAALKRKG LVQINNGERA RVS RPSADTI IGELSGMAKD FLAHPGGIAH FEQLRLFFES SLVRYAAENA TDAQIDLLAK ALEINSQSLD DNALFIRSDV DFHRVLAEIP GNPIFKAIHV ALLDWLIAAR PTVPDRELYE HNSVSYQQHI AIVDAIRQRD PDAADRALQT HLNSVSATWH ALGQQSKRKK
Source	Mammalian cell
Target Names	nanR
Protein Names	Recommended name: Transcriptional regulator NanR
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.