



Recombinant Microbacterium testaceum N-acyl homoserine lactonase (aiiM)

Product Code	CSB-MP512020MTN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C6L862
Product Type	Recombinant Protein
Immunogen Species	Microbacterium testaceum (strain StLB037)
Purity	≥85% (SDS-PAGE)
Sequence	MILAHDVSGS GPLLVLHGI TEDRRSWDPV DFTDGFTVVR VDLRGHGASA AEEPYDIPTL ATDVHDTLAQ LAENDVIPGE LPVIVGHSMG GIVATAYGAL FPARAIVNVD QPLQLAGMQG QVQQAEGMLR GADFPLFIHG MFAQMAGGLD AEELARVNGI RSPRQDVVLG MWRPLLEDSP EELAALVSGL TRIPEDVPYL VITGLDAGPE YAAWLQREIP QAVQEVWQPP THYPHLVDPA RFVERVEAFV R
Source	Mammalian cell
Target Names	aiiM
Protein Names	Recommended name: N-acyl homoserine lactonase Short name= AHL-lactonase EC= 3.1.1.81
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.