



Recombinant *Saccharomyces cerevisiae* Increasing suppression factor 1 (ISF1)

Product Code	CSB-EP512103SVN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C8ZET4
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MIASEIFERG VQDPFCQDCD YEDET DVQSF LGSNDL NDFV NSKLASF SFQ NSSKSNNSHH SSSTNAGNTS RHIGNHTIGH HLRKIKTAPH HLYGFVPANS TNNSNEPIRP SPRRIRANSS TLIHQLSRQS TRQSSLGDAA DSCFDHKCIK PRSRHSSCYG IPHLYGLEK YVSSELDLA VANDQSN DLT SPLTSVSTPA SNSNSYLNLN SSSAAYPSSY LSNEKNRRLK LISHGKISSN NVPGHSGNLN HYHRERTPSN LRRESFLLS NGSSSSPLQT RNNSYSNSLV KSPSNSSLNT SVASSNEESI PHTSNCLEER NPRRKSFIKL SLASSFSN
Source	E.coli
Target Names	ISF1
Protein Names	Recommended name: Increasing suppression factor 1 Alternative name(s): Mitochondrial biogenesis regulation protein 3
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.