



Recombinant *Saccharomyces cerevisiae* Restriction of telomere capping protein 4 (RTC4)

Product Code	CSB-EP512105SVN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C8ZFY3
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVGPGLGINR VRRKGVYSTK KGSGDNLLLM KRQGHDIHD RESDDLSGHD AFSPSKKRGK IDSITEDEIE VKKLSTVATF DKLSRSFPNS EVQAAKNAAL RGKEKEEEKV VSIPLIQNLK NEDIESIKCR NNLLLDGKKL LLEAELSAVE DNQIFSSSFP EDKKLSQLSC LSSKEQIIKK LQVREEYMSK FKLPPMLFSD ELLTEVEPFM PIVMDILEGK ISSAYYFEAK NAFKNSQKAY LSVDEFKRLN LNKFTAGFYG LKRQLRVGEE IAKRYKRALT HNQPATLKWV GITDFCNVYL APETLTSFCI YQLNLSNKSC SSKTPNKHPK QQLNEKEYYY DPELRMLAYD LLEDTVEYGI IVADSDPIEQ WEAAIEEDRL RELKLDVHNY SSRRWRLDTH D
Source	E.coli
Target Names	RTC4
Protein Names	Recommended name: Restriction of telomere capping protein 4
Expression Region	1-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.