



Recombinant *Saccharomyces cerevisiae* Enhancer of translation termination 1 (ETT1)

Product Code	CSB-YP513125SVL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C7GLH8
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain JAY291) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MAKRPLGLGK QSREKKRKVE SVEKKSDEPS RESTPVRSQM SVELDDADL DDELAQLKGL WSKYFHSDRD DEYVLNGIVH ECDRLRLSE EDKEIKKTLN DIFHGIYALA LSELTIFKAG DEEATEEKRK KDVSSFFESA IERVELGLSH FPESQFLKLV LAKIIFQRIP LEYISNLHLK SKDKKLDLVG QLEHGKHHFS IYENDTEFTF EILQMVNDLL DIVENFGREQ SIQEGIDSDN EEEEEIIDIE LEPEHPVYPL QQSLEANYEW LRNHFDKLLD NTNTDVKIYA SIANTLGELY LKKAEPSKV FLSLQYDDGG SEKVSDKEAK NAQETALKHT KKALEYLEKA KLEDDPDTWV QVAEAYIDLG NLLDNESAEQ EEAYKTAEEI LGKANKASHG KFQDVLDNFL QG
Source	Yeast
Target Names	ETT1
Protein Names	Recommended name: Enhancer of translation termination 1
Expression Region	1-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.