



Recombinant *Saccharomyces cerevisiae* DNA repair and recombination protein RDH54 (RDH54), partial

Product Code	CSB-BP513130SVL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C7GQI8
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain JAY291) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Source	Baculovirus
Target Names	RDH54
Protein Names	Recommended name: DNA repair and recombination protein RDH54 Alternative name(s): RAD homolog 54 Recombination factor TID1 Two hybrid interaction with DMC1 protein 1 Including the following 2 domains: DNA topoisomerase EC= 5
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.