



Recombinant *Saccharomyces cerevisiae* Altered inheritance of mitochondria protein 32 (AIM32)

Product Code	CSB-YP513133SVL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C7GS66
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain JAY291) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MLRITVKTLQ QRASFHHSFK HISVPDLHTR AQNDQTNCYC QEINARLPSK TDPLDPHIKL PHRTPNYNKH VLLSPGDRF AQPWKVAWNH NLDTNTNRPY NAISKLRSHL GGSPGILINA VHLQNEFIPR PKQHDEWLYF FVIPDMKLYV IKETDIEEFA SFLDEGAIQA PKLSFQDYLK GKAKASQQVH EVHHRKLTRF QGETFLRDWN LVCGHYKRDA KCGEMGPDII AAFQDEKLFP ENNLALISHI GGHIFAGNVI FYKLFGREKM QNKLDLWFG KVYPHNLKLL CENLENGKII DEMYRGGISM N
Source	Yeast
Target Names	AIM32
Protein Names	Recommended name: Altered inheritance of mitochondria protein 32
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.