



Recombinant *Saccharomyces cerevisiae* rRNA biogenesis protein RRP36 (RRP36)

Product Code	CSB-EP513138SVL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C7GWA5
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain JAY291) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSYYFKNLKP DLNSDVEEDD GNLLESIMAN KSKREIDEQE SSDDELKTLS FGSLKKAETI IDEEDFKDTK PVHKKPITTT YREESFDEDD DSEDKSEDA GFFEEDSEDE THHGQKVPKK KSKHAPIEQS SKKRVPVRN IPGLEIPRNL RSNLYQDIRF DKSTGKALDS SIIRKRYQFL DEYREKEIDE LQKLLQDRKF LSKIDQGERE EMEQRLKSMK SRLQSMKNKD LEREILKEYE NDMNKNNNTR YHLKKSEKRK VVQKWKFDHM KAKQREK VME RKRKRLGKE FKQFEFHNRR
Source	E.coli
Target Names	RRP36
Protein Names	Recommended name: rRNA biogenesis protein RRP36 Alternative name(s): Ribosomal RNA-processing protein 36
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.