



Recombinant *Saccharomyces cerevisiae* Biogenesis of lysosome-related organelles complex 1 subunit BLI1 (BLI1)

Product Code	CSB-YP513180SVN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C8ZCB1
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGEQNKLYYD VEKLVNSLQE SFDLDCAQSV SLFTSKRSRN EAWLEELNK FKLKDDVELD DVENLRAEID MKLNMLEDKV SYYERLYKEL EEFQNEIKIK TVVNNRRQSR TPK
Source	Yeast
Target Names	BLI1
Protein Names	Recommended name: Biogenesis of lysosome-related organelles complex 1 subunit BLI1 Short name= BLOC-1 subunit BLI1 Alternative name(s): BLOC-1 interactor 1
Expression Region	1-113
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.