



Recombinant Escherichia coli UPF0441 protein ygiB (ygiB)

Product Code	CSB-EP513292ENU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C4ZQV5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12 / MC4100 / BW2952)
Purity	>85% (SDS-PAGE)
Sequence	MKRTKSIRHA SFRKNWSARH LTPVALAVAT VFMLAGCEKS DETVSLYQNA DDCSAANPGK SAECTTAYNN ALKEAERTAP KYATREDCVA EFGEGQCQQA PAQAGMAPEN QAQAQQSSGS FWMPLMAGYM MGRLMGGGAG FAQQPLFSSK NPASPAYGKY TDATGKNYGA AQPGRMTMPV KTAMAPKPAT TTTVTRGGFG ESVAKQSTMQ RSATGTSSRS MGG
Source	E.coli
Target Names	ygiB
Protein Names	Recommended name: UPF0441 protein ygiB
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.