



Recombinant Escherichia coli Allantoinase (allB)

Product Code	CSB-EP513326ENU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C4ZUW1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12 / MC4100 / BW2952)
Purity	>85% (SDS-PAGE)
Sequence	MSFDLIKNG TVILENEARV VDIAVKGGKI AAIGQDLGDA KEVMDASGLV VSPGMVDAHT HISEPGRSHW EGYETGTRAA AKGGITTMIE MPLNQLPATV DRASIELKFD AAKGKLTIDA AQLGGLVSYN IDRLHELDEV GVVGFKCFVA TCGDRGIDND FRDVNDWQFF KGAQKLGELG QPVLVHCENA LICDELGEEA KREGRVTAHD YVASRPVFTE VEAIRRVLYL AKVAGCRLHV CHVSSPEGVE EVTRARQEGQ DVTCESEPHY FVLDTDQFEE IGTLAKCSPP IRDLENQKGM WEKLFNGEID CLVSDHSPCP PEMKAGNIMK AWGGIAGLQS CMDVMFDEAV QKRGMSLPMF GKLMATNAAD IFGLQQKGRI APGKDADFVF IQPNSSYVLT NDDLEYRHKV SPYVVRTIGA RITKTILRGD VIYDIEQGFP VAPKGQFILK HQQ
Source	E.coli
Target Names	allB
Protein Names	Recommended name: Allantoinase EC= 3.5.2.5 Alternative name(s): Allantoin-utilizing enzyme
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.