



Recombinant Escherichia coli Isocitrate dehydrogenase kinase/phosphatase (aceK)

Product Code	CSB-YP513369ENU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C5A0V3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12 / MC4100 / BW2952)
Purity	>85% (SDS-PAGE)
Sequence	MPRGLELLIA QTILQGFDAQ YGRFLEVTSQ AQRFEQADW HAVQQAMKNR IHLYDHHVGL VVEQLRCITN GQSTDAFLL RVKEHYTRLL PDYPRFEIAE SFFNSVYCRL FDHRSALTER LFIFSSQPER RFRTIPRPLA KDFHPDHGWE SLLMRVISDL PLRLRWQNKSRDIHYIIRHL TETLGTDNLA ESHLQVANEL FYRNKAAWLVGKLTIPSGTL PFLPIHQTD DGELFIDTCL TTTAEASIVF GFARSYFMVY APLPAALVEW LREILPGKTT AELYMAIGCQ KHAKTESYRE YLVLQGCNE QFIEAPGIRG MVMLVFTLPG FDRVFKVIKD RFAPQKEMSA AHVRACYQLV KEHDRVGRMA DTQEFENFVL EKRHISPALM ELLLQEAAEK ITDLGEQIVI RHLYIERRMV PLNIWLEQVE GQQLRDAIEE YGNAIRQLAA ANIFPGDMLF KNFGVTRHGR VVFYDYDEIC YMTEVNFRDI PPRYPEDDEL ASEPWYSVSP GDVFPEEFRH WLCADPRIGP LFEEMHADLF RADYWRALQN RIREGHVEDV YAYRRRQRFS VRYGEMLF
Source	Yeast
Target Names	aceK
Protein Names	Recommended name: Isocitrate dehydrogenase kinase/phosphatase Short name= IDH kinase/phosphatase Short name= IDHK/P EC= 2.7.11.5 EC= 3.1.3.-
Expression Region	1-578
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.