



# Recombinant *Saccharomyces cerevisiae* Altered inheritance of mitochondria protein 32 (AIM32)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP513911SVN   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | C8ZEF8  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MLRITVKTLQ QRASFHHSFK HISVPDLHTR AQNDQTNCYC QEINARLPSK<br>TDPLDPHIKL PHRTPNYNKH VLLLSPGDRF AQPWKVAWNH NLDNTNRPY<br>NAISKLRSHL GGSPGILINA VHLQNEFIPR PKQHDEWLYF FVIPDMKLYV<br>IKETDIEEFA SFLDEGAIQA PKLSFQDYLS GKAKASQQVH EVHHRKLTRF<br>QGETFLRDWN LVCGHYKRDA KCGEMGPDII AAFQDEKLPF ENNLALISHI<br>GGHIFAGNVI FYKLFGREKM QNKLDLWFG KVYPHNLKLL CENLENGKII<br>DEMYRGGISM N                                      |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | AIM32   |
| <b>Protein Names</b>     | Recommended name: Altered inheritance of mitochondria protein 32  |
| <b>Expression Region</b> | 1-311   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |