



Recombinant Escherichia coli Tetraacyldisaccharide 4'-kinase (lpxK)

Product Code	CSB-EP514353ENU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C4ZQ41
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12 / MC4100 / BW2952)
Purity	>85% (SDS-PAGE)
Sequence	MIEKIWSGES PLWRLLLPLS WLYGLVSGAI RLCYKCLKLR AWRAPVPVVV VGNLTAGGNG KTPVVVWLVE QLQQRGIRVG VVSRGYGGKA ESYPLLLSAD TTTAQAGDEP VLIYQRTDAP VAVSPVRSDA VKAILAQHPD VQIIVTDDGL QHYRLARNVE IVVIDGVRRF GNGWWLPAGP MRERAGRLKS VDAVIVNGGV PRSGEIPMHL LPGQAVNLRG GTRCDVAQLE HVVAMAGIGH PPRFFATLKM CGVQPEKCVL LADHQLSLNHA DVSALVSAGQ TLMVTEKDAV KCRAFAEENW WYLPVDAQLS GDEPAKLLTQ LTLLASGN
Source	E.coli
Target Names	lpxK
Protein Names	Recommended name: Tetraacyldisaccharide 4'-kinase EC= 2.7.1.130 Alternative name(s): Lipid A 4'-kinase
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.