



# Recombinant Escherichia coli Anhydro-N-acetylmuramic acid kinase (anmK)

<b>Product Code</b>	CSB-YP514408ENU
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	C4ZYA3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12 / MC4100 / BW2952)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKSGRFIGVM SGTSLDGV DV VLATIDEHRV AQLASLSWPI PVSLKQAVLD ICQQQQLTLS QFGQLDTQLG QLFADAVNAL LKEQNLQARD IVAIGCHGQT VWHEPTGVAP HTLQIGDNNQ IVARTGITVV GDFRRRDIAL GGQGAPLVPA FHHALLAHPT ERRMVLNIGG IANLSLLIPG QPVG GYDTGP GNMLMDAWIW RQAGKPYDKD AEWARAGKVI LPLLQNMLSD PYFSQPAPKS TGREYFN YGW LERHLRHFP G VDPRDVQATL AELTAVTISE QVLLSGGCER LMVCGGGSRN PLLMARLAAL LPGTEVTTTD AVGISGDDME ALAFAWLAWR TLAGLPGNLP SVTGASQETV LGAIFPANP
<b>Source</b>	Yeast
<b>Target Names</b>	anmK
<b>Protein Names</b>	Recommended name: Anhydro-N-acetylmuramic acid kinase EC= 2.7.1.170 Alternative name(s): AnhMurNAc kinase
<b>Expression Region</b>	1-369
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.