



Recombinant Emericella nidulans Autophagy-related protein 3 (atg3)

Product Code	CSB-MP514955EKF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C8VDI2
Product Type	Recombinant Protein
Immunogen Species	Emericella nidulans (strain FGSC A4 / ATCC 38163 / CBS 112.46 / NRRL 194 / M139) (Aspergillus nidulans)
Purity	>85% (SDS-PAGE)
Sequence	MNILHSTLST WRDRLAPVSR TSTFRTTGQI TPEEFVLAGD YLVYKFPSWS WGDASSPSKR VSYLPPGKQF LVTRGVPCHR RLNENFAGDA HLDDEIVRDF LSGGAGDSDG GDDNDGWLRT GGGGKRHEST IRDVRTVDES GKEEAEEVEE EDIPDMEDDD DEEAIREPA GTTSTTQPTR TYNLYITYSN FYRTPRLYLS GYLSPSEPLP PKLMMEDIVG DYKDKTVTLE DFPWFDGSLQ MASVHPCRHA SVMKTLLDRA DAALKLRREK IKQTAASSSP QEKAKLAQPE SGLEGLVDDI KGLSLGDAQK QCQGDKGQAG GDEWEVLQHD EEEVAIRVD QYLVVFLKFI ASVTPGIEHD FTMGV
Source	Mammalian cell
Target Names	atg3
Protein Names	Recommended name: Autophagy-related protein 3 Alternative name(s): Autophagy-related E2-like conjugation enzyme atg3
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.