



# Recombinant *Saccharomyces cerevisiae* KRR1 small subunit processome component (KRR1)

<b>Product Code</b>	CSB-MP514959SVN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	C8Z430
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVSTHNRDKP WDTDDIDKWK IEEFKEEDNA SGQPFAEESS FMTLFPKYRE SYLKTIWNDV TRALDKHNIA CVL DLVEGSM TVKTTRKTYD PAILKARDL IKLLARSVPF PQAVKILQDD MACDVIKIGN FVTNKERFVK RRQRLVGPNG NTLKALELLT KCYILVQGNT VSAMGPFKGL KEVRRVVEDC MKNIHPIYHI KELMIKRELA KRPELANEDW SRFLPMFKKR NVARKPKKI RNVEKKVYTP FPPAQLPRKV DLEIESGEYF LSKREKQMKK LNEQKEKQME REIERQEERA KDFIAPEEEA YKPNQN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KRR1
<b>Protein Names</b>	Recommended name: KRR1 small subunit processome component Alternative name(s): KRR-R motif-containing protein 1 Ribosomal RNA assembly protein KRR1
<b>Expression Region</b>	1-316
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.