



Recombinant *Saccharomyces cerevisiae* Spindle pole component BBP1 (BBP1)

Product Code	CSB-EP515236SVK-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	E7NNK3
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain FostersO) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNQEDNTGGG GIFGLFKWTK XALFGTDISP SMKYKDQEER RDRSRYAQDD TNFSMKFGND SNRRSTNLSR SNSWSGLDST LHRKYELLPE YNENGFNSIV NGDHHSKERI RSLRSPAPIV PREPLRNEPT DTFGHLHTK RRTINELSNS QIPFIPPQED DPLLSKLFNK DGVNEVRRSP YKLSVKDIPG KFPSPLTKRD EIDNYYVRDE DACHKNREYK KAYFDLFAQM DLNSRDLEDL CEDVREQREQ FHRNEQTYKQ AYEEMRAELV NELKSKTLF ENYYSLGQKY KSLKKXLDQT ISHEAELATS RERLYQEEDL KNFEIQTLKQ RLSDLELKYT NLQIEKDMQR DNYESEIHDL LLQLSLRNNE RKDTSAGSNI FSTGQYDRTP FHNGNNSYDS NSHSWDTDYL KKYRRIHRTL KQKS
Source	E.coli
Target Names	BBP1
Protein Names	Recommended name: Spindle pole component BBP1 Alternative name(s): BFR1-binding protein 1
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.