



# Recombinant Tropomodulin (unc-94)

<b>Product Code</b>	CSB-MP515310CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O01479
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSQAKTDYYS EEKTFAPS SA NSQQGTQLPS KVYNKGLKDL EDNDIEGLLS SLSIDELEDL NNDFDPDNSM LPPSQRCDQ TDKEPTGPYK RDNLLKFLED KAKTEKDVED VCPYTPGQKR GKVYDSDSGR NSEEPENKGM EMPIEIDLDD DEEELECALV TAPEKDLVDL AGILGMHNVL NQPQYYNALK GKTQDESTGT TFNGIMQSYV PRIVPDEPDN DTDVESCINR LREDDTDLKE VNINNMKRVS KERIRSLIEA ACNSKHIEKF SLANTAISDS EARGLIELIE TSPSLRVLVN ESNFLTPELL ARLLRSTLVT QSIVEFKADN QRQSVLGNQV EMDMMAIEE NESLLRVGIS FASMEARHRV SEALERNYER VLRRLGKDP NV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	unc-94
<b>Protein Names</b>	Recommended name: Tropomodulin Alternative name(s): Tmod-like protein Short name= cTmd1 Uncoordinated protein 94
<b>Expression Region</b>	1-392
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.