



Recombinant Multiple exostoses homolog 1 (rib-1)

Product Code	CSB-BP515314CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O01704
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MKFHLVIFML FGSVRLQNPT IERKQCTMSN CFDFSKCSTS KKVYIHPMEK RFEESPQSVI YSKILKHFL E SNHYTNDPNE ACIFLLGIDT TDRDVRSQNY VKNVNDYIES LDPSVWNNGR NHLIFNFYHG TFPDYDDHNL NFDTGEMIA RASSENFI KVFDVSLPLF HENHPYEIKE SKSERNDDRI ENQRKYLVSF K GKRYVYGIG SGTRNLVHHL HNGDDIVMVT TCKHNNDWQV YQDDRCQRDN DEYDRWEYDE LLANSTFCLV PRGRRLGSFR FLETLRSGCV PVVISDSWIL PFSETIDWNS AAIVVAERDA LSIPELLMST SRRRVKELRE SARNVYDAYL RSIQVISDHV LRIIFKRIDN KIELEDHQ
Source	Baculovirus
Target Names	rib-1
Protein Names	Recommended name: Multiple exostoses homolog 1 Alternative name(s): Related to mammalian RIB protein 1
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.