



Recombinant *Schizosaccharomyces pombe* Uncharacterized WD repeat-containing protein C25H1.08c (SPAC25H1.08c)

Product Code	CSB-MP515617SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13982
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MAEEEEENQEL YIAENEVEEV IKDNDK LMDQ GVTEDEAEQE QNDTMN NMVM DMSVQGGFFEH KDSVFSVSIN PVHSNLCASG GGDDLGYIWD ITTGEQICQL TGHKDSIVAI DWSFDGTYIA TGGMDSQVRL WKSSTGFEFI TAFETVDEIV WLSWHPKGLF LAAGCNDGSI WMWSLPSGKV VQVMYGHTAP VNAGKFIPPG VGKRLATVDD SGT LIVWNP A TGAPECRMSS DDHRFDPGNE ETAAGWTSFD CNAEGNVLFL GGSSGKVKVV NINSSHILAS LEAQTESVEA IALCTALPIC ACASVDGTVA LYDSASLKFR KSLPHEQAVI DCKFLPNTPY LLTACADCVI RKWDVRSQQL LGEYTG HQEP ILCMAITPDG KRVVTGSDDT ELLVFDCEH
Source	Mammalian cell
Target Names	SPAC25H1.08c
Protein Names	Recommended name: Uncharacterized WD repeat-containing protein C25H1.08c
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.