



Recombinant Schizosaccharomyces pombe Pre-mRNA-splicing factor cwc21 (cwf21)

Product Code	CSB-YP515646SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14161
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSYNGIGLPT PRGSGTNGYV MRNLSHVKKY DKNTNLQSNR NAKALEKRVQ DPSISEHECR RQIESKLLLY REQLLEEVS QHSTDAAASD SNTNFGTENP KPPKAIKDE QSQSKTKSLD EADVEILVQK YREQLLKEQ LQKSTEKGN FESILQPQKR KETRGFHSKN DDDGRLEERD DLRSSVDDDY YDYPRYSERK SLNSKRHVDN YNENRRRHYD SYSSYDELER RRSSNESYSR RSELPRRDYN RHDERERYSY HRRRERSNSP SYTKNESIPV VDRDSSPEGG EIV
Source	Yeast
Target Names	cwf21
Protein Names	Recommended name: Pre-mRNA-splicing factor cwc21 Alternative name(s): Cell cycle control protein cwf21 Complexed with cdc5 protein 21
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.