



Recombinant Mouse Galectin-3 (Lgals3)

Product Code	CSB-MP012887MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P16110
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	ADSFSLNDA LAGSGNPNPQ GYPGAWGNQP GAGGYPGAAY PGAYPGQAPP GAYPGQAPPG AYPGQAPPSA YPGPTAPGAY PGPTAPGAYP GQPAPGAFPG QPGAPGAYPQ CSGGYPAAGP YGVPAAGPLTV PYDLPLPGGV MPRMLITIMG TVKPNANRIV LDFRRGNDVA FHFNPRFNEN NRRVIVCNTK QDNNWGKEER QSAFPFESGK PFKIQVLVEA DHFQVAVNDA HLLQYNHRMK NLREISQLGI SGDITLTSAN HAMI
Source	Mammalian cell
Target Names	Lgals3
Protein Names	Recommended name: Galectin-3 Short name= Gal-3 Alternative name(s): 35 kDa lectin Carbohydrate-binding protein 35 Short name= CBP 35 Galactose-specific lectin 3 IgE-binding protein L-34 galactoside-binding lectin
Expression Region	2-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.