



Recombinant Schizosaccharomyces pombe SWR1-complex protein 4 (swc4)

Product Code	CSB-MP515670SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14308
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MTSADIRDVF ELPPPEIGNK QKSKTPTErr PEGISRELYS LLGENSAPLA IYQKKFKEKP KVSHKAKNWV RQPFSSSRK DDFTLHHWVL KSEVDSEASY KFEKFNVPLF IIDYTDEEYQ NYLKDEDWNK DETDYLFRLC KDYDLRFFVI ADRYDNEKYK KHRTLEDLKD RFYSVSRKIL LARNPINSMT AAQSSLLNTM EYNKEQEVIR KKYLIGLASR TPEEVAEEEE LFIELKRIET SQAKLLSDRD EVLRLLEQK GDGGIHEYHT SAGMSSLIQD MINSQRTKNK VEEAIVSSSA PSSGVSSVLN TPTRPHALST PRIRYGPQPT DPQFGITWHE KLHPGTFVRS QKIPAİKASL SQRVSSVMTE LGVSSRLIMP TAKNFEKFVE LQNSIVSLLI LKRKVDRLSQ ETEIQDKLSR KRSASPDGSE SKKHISQ
Source	Mammalian cell
Target Names	swc4
Protein Names	Recommended name: SWR1-complex protein 4
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.