



# Recombinant *Schizosaccharomyces pombe* Serine/threonine-protein kinase ksp1 (ksp1)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-BP515673SXV  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | O14328   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | <i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | MKLLQKKGYK VERPLNKGSY GTVVLAHRLF RTPRCKDLKY AIKCIKKPAY<br>TFLQEVNILR QLSRSRHRNI IHVESFEDN VYYYVVLEYC PLGDLYECIL<br>NNDFPNAKNQ PEMIKNIFLQ IIDGVAHLHS HGIYHRDLKP ENFLLSLED<br>GSELVVKISD FGLACRDKIS YDFGTGSDRY MAPEQFEEVD GAGYSPRAAD<br>IWALGICLLN LIFARNPPTY PHEKDPIFAD YMLDAMTLFD VFPTLSQDTY<br>NVLRACLCVS PEKRSKAKTR EAVLAVTKWT TDDEELESFV NEEEEFRASD<br>FMPAEDNVRC TQSDREPLRT PSVLTPANTI QRGLLPSKLP ALSDV DENIS<br>TSSSPRSPAS LAPVNNERS YDSGLGESLN NMHIGKSIAT AVPVNTRKSP<br>YSCSAPAIVF PNSIKGNKDH LKFGRSWCDM DEEDEEDIVS FGSNDDFGAS<br>DELSSKHIGL ADDWNVLSQW NDNS |
| <b>Source</b>            | Baculovirus  |
| <b>Target Names</b>      | ksp1   |
| <b>Protein Names</b>     | Recommended name: Serine/threonine-protein kinase ksp1 EC= 2.7.11.1<br>Alternative name(s): Serine/threonine-protein kinase ppk20  |
| <b>Expression Region</b> | 1-474  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |