



Recombinant *Drosophila melanogaster* Basic helix-loop-helix neural transcription factor TAP (tap)

Product Code	CSB-MP515701DLU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O16867
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAACYNAYSA GSQSFEFDED DDDASFDSGY EKSFETEAQL SSRRRLDFGT PPTPAIPQPY SGGTWDVAVPL SSPPAGFVGL LDTSSNHSTR SGRTLVEHLN SRATNGVFDP PLTSTPVKSP EDPNAPRPKR KYAVGKNRVT RSRSPQVVK IKRFRMKAN DRERNRMHNL NDALEKLRVT LPSLPEETKL TKIEILRFAH NYIFALEQVL ESGGSINLDL EKLQNFTLSG ERITKELFDA LFNVPQPYPL FGRMFYGGG MAPLAQHQA PASHAEQPPA MGGFQHGMDY PQQPPGFDFD GSMRFYHQQQ QQPHQPHHLQ PNPQQESSPQ QFSQEKYDLF RGSFDAAANL HSTNLDSGIH QQSSFYSQTP PWKDYPEDQA HVHPVPHQHS YKNFAPQV
Source	Mammalian cell
Target Names	tap
Protein Names	Recommended name: Basic helix-loop-helix neural transcription factor TAP Alternative name(s): Protein biparous Target of Poxn protein
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.