



Recombinant *Drosophila ananassae* Alpha-amylase 2 (Amy58)

| | |
|--------------------------|---|
| Product Code | CSB-EP515721DKY |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | O18345 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila ananassae</i> (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Sequence | <p>QF NTNYASGRSG MVHLFEWKWD DIAAECENFL GPYGYAGVQV SPVNAVND SRPWWERYQP ISYKLVTRSG NEEQFASMVR RCNNVGVRIY VDVVFNHMAA DGGTYGTGGS TASPSSKSYG GVPFSSLDNF PTCAISNYND ANQVRNCELV GLRDLNQQNS YVQEKIVEFL NHLIDLGVAG FRVDAAKHMW PADLGVIYGS LKNLNTDHGF ESGAKAYIVQ EVIDMGGEAI SKSEYTLGA ITEFRHSDSI GKAFRGKNQL QYLVNWGVSF GFAASDRSLV FVDNHDNQRG HGAGGADVLT YKVPKQYKMA SAFMLAHPFG TPRVMSSFSF TDTDQGPPTT DGQNIASPSF NSDNSCSGGW VCEHRWKQIY NMVGFRNAVQ SDAIQNWWDN GSNQIAFSRG SKGFVAFNND NYDLNSSVQT GLPAGTYCDV ISGSKSGSSC TGKTVTVGSD GRANISIGSS EDDGVLAIHV NAKL</p> |
| Source | E.coli |
| Target Names | Amy58 |
| Protein Names | Recommended name: Alpha-amylase 2 EC= 3.2.1.1 |
| Expression Region | 19-494 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |