



# Recombinant Arabidopsis thaliana GDSL esterase/lipase At2g30310 (At2g30310)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP515794DOA-B   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | O22927  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Arabidopsis thaliana (Mouse-ear cress)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | QP LFPAILIFGD STVDTGNNNY HSQTIFKAKH LPYGVLDLPGH<br>EANGRYSNGK VISDVIASKL NIKELVPPFL QPNISHQDIV TGVSFASAGA<br>GYDDRSSLSS KAIPVSQQPS MFKNYIARLK GIVGDKKAME IINNALVVIS<br>AGPNDFILNF YDIPTRRLEY PTIHGYQEFI LKRLDGFVRE LYSLGCRNIV<br>VGGLPPMGCL PIQMTAKMRN ILRFCVEQEN KDSVLYNQKL VKKLPEIQAS<br>LPGSNFLYAN VYDPLMDMIQ NPSKYGFKET KKGCCGTGYL ETTFMCNPLT<br>KTCPNHSDHL FWDSIHPSEA AYNYIGNFVD AQIRGWIKA             |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | At2g30310   |
| <b>Protein Names</b>     | Recommended name: GDSL esterase/lipase At2g30310 EC= 3.1.1.- Alternative name(s): Extracellular lipase At2g30310  |
| <b>Expression Region</b> | 29-359  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |