



Recombinant Solanum tuberosum Cysteine protease inhibitor 8

Product Code	CSB-MP515835FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O24384
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	LVLPEVYD QDGNPLRIGE RYIIKNPLLG AGAVYLDNIG NLQCPNAV LQ HMSIPQFLGK GTPVVFIRKS ESDYGDVVRL MTAVYIKFFV KTTKLCVDET VWKVNNEQLV VTGGNVGNEN DIFKIKKTDL VIRGMKNVYK LLHCPSHLEC KNIGSTFKNG YPRLVTVNDE KDFIPFVFIK A
Source	Mammalian cell
Protein Names	Recommended name: Cysteine protease inhibitor 8 Alternative name(s): PCPI-8 Short name= PcpI8
Expression Region	43-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.