



Recombinant Human Histo-blood group ABO system transferase (ABO), partial

Product Code	CSB-EP001110HU1-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P16442
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AVREPDHLQRVSLPRMVYPQPKVLTPCRKDVLVVTWLPWLAPIVWEGTFNIDILN EQFRLQNTTIGLTVFAIKKYVAFLKLFLETAEKHFVMVGHVHYYVFTDQPAAVP RVTLTGTGRQLSVLEVRAYKRWQDVSMRRMEMISDFCERRFLSEVDYLVCVDV DMEFRDHVGV EILTPLFGTLHPGFYGS SREAFTYERRPQS QAYIPKDEGDFY LGGFFGGSVQEVQRLTRACHQAMMVDQANGIEAVWHDESHLNKYLLRHKPT KVLSP EYLWDQQLLGWPAVLRKLRFTAVPKNHQAVRNP
Source	E.coli
Target Names	ABO
Protein Names	Recommended name: Histo-blood group ABO system transferase Alternative name(s): Fucosylglycoprotein 3-alpha-galactosyltransferase Fucosylglycoprotein alpha-N-acetylgalactosaminyltransferase Glycoprotein-fucosylgalactoside alpha-N-acetyl gal
Expression Region	54-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This gene encodes proteins related to the first discovered blood group system, ABO. Which allele is present in an individual determines the blood group. The O blood group is caused by a deletion of guanine-258 near the N-terminus of the protein which results in a frameshift and translation of an almost entirely different protein. Individuals with the A, B, and AB alleles express glycosyltransferase activities that convert the H antigen into the A or B antigen. Other minor alleles have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.