



Recombinant *Shigella flexneri* Outer membrane protease icsP (icsP)

Product Code	CSB-BP516443SZB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O33641
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	<i>Shigella flexneri</i>
Purity	>85% (SDS-PAGE)
Sequence	TTNYPLFIPD NISTDISLGS LSGKTKERVY HPKEGGRKIS QLDWKYSNAT IVRGGIDWKL IPKVSFGVSG WTTLGNQKAS MVDKDWNNNSN TPQVWTDQSW HPNTHLRDAN EFELNLKGWL LNNLDYRLGL IAGYQESRYS FNAMGGSYIY SENGGSRNKK GAHPSGERTI GYKQLFKIPY IGLTANYRHE NFEFGAELKY SGWVLSSDTD KHYQTETIFK DEIKNQNYCS VAANIGYYVT PSAKFYIEGS RNYISNKKGD TSLYEQSTNI SGTIKNSASI EYIGFLTSAG IKYIF
Source	Baculovirus
Target Names	icsP
Protein Names	Recommended name: Outer membrane protease icsP EC= 3.4.23.-
Expression Region	21-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.