



Recombinant Schizosaccharomyces pombe Anaphase-promoting complex subunit 14 (apc14)

Product Code	CSB-EP516644SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42659
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDPFAGTGQS KSYRTFGNVF QRLSMSGNPK PLAKKPLLLP EASKAIEFWN YRDVIGGEEA EQERNERASR IAISRIASSR SSIPEYRQAI MQHLTKHVQQ LDQLAEW
Source	E.coli
Target Names	apc14
Protein Names	Recommended name: Anaphase-promoting complex subunit 14 Alternative name(s): 20S cyclosome/APC complex protein apc14 Overlapping meiotic transcript protein 1
Expression Region	1-107
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.