



Recombinant FKBP-type peptidyl-prolyl cis-trans isomerase fkpA (fkpA)

Product Code	CSB-BP517042AXW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O08437
Product Type	Recombinant Protein
Immunogen Species	Aeromonas hydrophila
Purity	>85% (SDS-PAGE)
Sequence	CQKDEKTAAN TAEVKAEASK PAEAPKAEAK SFEEQSGYAI GLSMGRYIAN TLERQQELGI KLDNSVILKG VTDGLGKEAK MTDEEIQKVL QQYDAKINEL TKAKADKDAV ENQKKGEEYL AANAKKEGVK STESGLQYQV EKMGTGAKPK ATDIVKVHYT GTLTDGTKFD SSVDRGEPAT FPLNQVIPGW TEGVQLMPVG SKFKFFLPSK LAYGEHGAGS IPANAVLVFD VELLAIEKPA ADGDNAKK
Source	Baculovirus
Target Names	fkpA
Protein Names	Recommended name: FKBP-type peptidyl-prolyl cis-trans isomerase fkpA Short name= PPIase EC= 5.2.1.8 Alternative name(s): Rotamase
Expression Region	21-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.