



# Recombinant *Saccharomyces cerevisiae* Topoisomerase I damage affected protein 11 (TDA11)

<b>Product Code</b>	CSB-EP517215SVK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	E7NIP9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain FostersO) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNKFDEFIES NEKDLVDVTS TRNSIISMSP VRKTGRKIRS ASSNGYRLEH HRTSSAGSMH SQRLMTPTRL NDQDHPLQAK PDARRVVTRH SSVSVPNAMS KRRSLIQPMV VPTTPESQNN LPFGGSVSHS EGSYGIPLES TTVLSSEQAM ASGLRRSRNG SSQSVNSMIA TTIPTNGVDV SALLQSLATK ELELLECKQK IEDLKKQTQH EEQNYTRRAC ELHELKEQVS KHLDP SLNTP VKNRAFSPVY QNIPLES RTE NAGNSSLPSS VSKPKNMGHQ STNQSR SVSP QDIQERRQRD DSSDSSKQSL WSKPLALFNQ FDKIIQHEIE RTLNWDD SLS GTPEVQEGTP TSNSESSAQQ YDNEAPGARQ KSPSQGSVSR SLWSFVSDVK AGLLGIEEEN DNDVITDNRC DPVYKSDRQH EQKKSTHKIT NRGQAEDSGD DSSLNMRKFK TTKFQKDNA GNNSLTDESG HRTREKKK
<b>Source</b>	E.coli
<b>Target Names</b>	TDA11
<b>Protein Names</b>	Recommended name: Topoisomerase I damage affected protein 11
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.