



Recombinant *Saccharomyces cerevisiae* Mitochondrial morphogenesis protein SLD7 (SLD7)

Product Code	CSB-YP517221SVR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	E7QKN3
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Zymaflore VL3) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSRKLCTLNF TLSGKQGSLV IRDIQLWSNR PTASKSTSEL RGQFIQYVDL AKLPLWVRST NMNTYRCYST SATAQAYFKS KLRNANRGIV IELSDKVDQR SQEPAYLIIF RENTELNCFQ VDLTMKHEFD GQVTKLKQEI GKTRASVSKE GSIDIIQQS QQRKIGTKTK VYRNVHINDK RLQFNETLSK LILGGLRLRG ISNSITDYQK LYKITFDAAE FTHRDELKRI SMGSGEEVSF ESLQETVETL LKLFTKS
Source	Yeast
Target Names	SLD7
Protein Names	Recommended name: Mitochondrial morphogenesis protein SLD7 Alternative name(s): Synthetic lethality with DPB11-24 mutation protein 7
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.