



# Recombinant Arabidopsis thaliana Ribulose bisphosphate carboxylase large chain (rbcL)

<b>Product Code</b>	CSB-YP517305DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O03042
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	PQTETKAS VGFKAGVKEY KLTYTPEYE TKDTDILAAF RVTPQPGVPP EEAGAAVAE SSTGTWTTVW TDGLTSLDRY KGRCYHIEPV PGEETQFIAY VAYPLDLFEE GSVTNMFTSI VGNVFGFKAL AALRLEDLRI PPAYTKTFQG PPHGIQVERD KLNKYGRPLL GCTIKPKLGL SAKNYGRAVY ECLRGGLDFT KDDENVNSQP FMRWRDRFLF CAEAIYKSQA ETGEIKGHYL NATAGTCEEM IKRAVFAREL GVPIVMHDYL TGGFTANTSL SHYCRDNGLL LHIHRAMHAV IDRQKNHGMH FRVLAKALRL SGGDHIHAGT VVGKLEGDRE STLGFVDLLR DDYVEKDRSR GIFFTQDWVS LPGVLPVASG GIHVWHMPAL TEIFGDDSVL QFGGGTLGHP WGNAPGAVAN RVALEACVQA RNEGRDLAVE GNEIIREACK WSPELAAACE VWKEITFNFP TIDKLDGQE
<b>Source</b>	Yeast
<b>Target Names</b>	rbcL
<b>Protein Names</b>	Recommended name: Ribulose bisphosphate carboxylase large chain Short name= RuBisCO large subunit EC= 4.1.1.39
<b>Expression Region</b>	3-479
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.