



Recombinant Cynodon dactylon Profilin (PRO1)

Product Code	CSB-MP517332EQB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O04725
Product Type	Recombinant Protein
Immunogen Species	Cynodon dactylon (Bermuda grass) (Panicum dactylon)
Purity	≥85% (SDS-PAGE)
Sequence	SWQAYVDDH LMCEIEGHHL TSAAIIGHDG TVWAQSAAFP AFKPEEMANI MKDFDEPGFL APTGLFLGPT KYMVIQGEPE AVIRGKKKSG GVTVKKTGQA LVIGIYDEPM TPGQCNMVIE KLG DY LIEQG M
Source	Mammalian cell
Target Names	PRO1;
Protein Names	Recommended name: Profilin Alternative name(s): Pollen allergen Cyn d 12 Allergen= Cyn d 12
Expression Region	2-131
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.