



Recombinant Arabidopsis thaliana Probable histone H2AXa (At1g08880)

Product Code	CSB-EP517333DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O04848
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSTGAGSGTT KGGRGKPKAT KSVSRSSKAG LQFPVGRIAR FLKSGKYAER VGAGAPVYLS AVLEYLAAEV LELAGNAARD NKKTRIVPRH IQLAVRNDEE LSKLLGSVTI ANGGVLPNIH QTLPSKVGK NKGDIGSASQ EF
Source	E.coli
Target Names	At1g08880
Protein Names	Recommended name: Probable histone H2AXa Alternative name(s): HTA5
Expression Region	1-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.