



Recombinant *Thermotoga maritima* Cell division protein FtsZ (ftsZ)

Product Code	CSB-MP517482TNJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O08398
Product Type	Recombinant Protein
Immunogen Species	<i>Thermotoga maritima</i> (strain ATCC 43589 / MSB8 / DSM 3109 / JCM 10099)
Purity	>85% (SDS-PAGE)
Sequence	MGFDLDVEKK KENRNIPQAN NLKIKVIGVG GAGNNAINRM IEIGIHGVEF VAVNTDLQVL EASNADVKIQ IGENITRGLG AGGRPEIGE QAALESEEKIR EVLQDTHMVF ITAGFGGGTG TGASPVIKI AKEMGILTVA IVTTPFYFEG PERLKAIEG LKKLRKHVDT LIKISNNKLM EELPRDVKIK DAFLKADETL HQQVKGISEL ITRGYINLD FADIESVMKD AGAAILGIGV GKGEHRAREA AKKAMESKLI EHPVENASSI VFNITAPSNI RMEEVHEAAM IIRQNSSEDA DVKFGLIFDD EVPDDEIRVI FIATRFPEDED KILFPEGDIP AIYRYGLEGL L
Source	Mammalian cell
Target Names	ftsZ
Protein Names	Recommended name: Cell division protein FtsZ
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.