



Recombinant Kluyveromyces lactis Protein SPT3 (SPT3)

Product Code	CSB-MP517535K BK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13472
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MSDKYKYRVE IQQMMFVSGE NNEPPVETTS LIEDIVRGQV VEILLQATRT ANCRGSKSIV PEDVIFLIRH DKAKVNRLRT YLSWKDLRKN AKDQDAQQAA AGGGAGDGS G PGGIDDDAG KKDKDGAGVN MKVKKSSIKL PWELQFMFNE QPLENPDEDD MDDDEREANM ATLKRLKMAD DRTRGMTKEE YVHWSDCRQA SFTFRKSKRF KDWSGISQLI DSKPHDDVID ILGFLTFEIV CAITETALKI KTRYERIRSM KTQSQPQEYS VGTTRRKRL FDGPDNVVNP LKPEHIEEAW RVLQTVNMKH RALSNYKGGR LSSRTVIM
Source	Mammalian cell
Target Names	SPT3
Protein Names	Recommended name: Protein SPT3
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.