



Recombinant Schizosaccharomyces pombe Protein pof8 (pof8)

Product Code	CSB-EP517575SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13795
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MFVPRQLNVR KIKAFGTGKEN NSIADGNNNK LKDEHYKHNE ASKEPSHSIS GGLMLNQQDR QLIEPLNPDF LSAVDSILEI YFHRERQKEK VHFLAFLIQQD DFWKGIRPNP TQNNLKYALS YVTNALFHFD NSSHMVIRNE NIVLPLDIPL YDRIIYVEPV PATLSNKSL LAGKLRKYLK EFLPYVDAIG TPEGYAFVIL YKKVDQSALS KLPVPPGWVL LTRKEWTNRE EKYFENQLHL VKASSSDVSN SSNSFPENRY PKLTKVEKQM TKSVSQTSQT DKDEDNLDFT KNLTRIKNL HPLTNKSTIH SLLSYVFSRQ TQNIACEPMY IDYRKDETEA IIRWKTP LHA ETCINAFRTQ ERKQNSHDDI RAHRKKGSSR PFLIAELITG EEEKNYWRML KK
Source	E.coli
Target Names	pof8
Protein Names	Recommended name: Protein pof8
Expression Region	1-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.