



Recombinant Schizosaccharomyces pombe U1 small nuclear ribonucleoprotein 70 kDa homolog (usp101)

Product Code	CSB-EP517581SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13829
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAEKLPAPLL ALFAPRPPLR YLPPMDVPPE KRSTPRVSGI AKYLKYAQSH DQQYHPTESL EEKRLRLRDE KQKQQRERLR SMIKVVDPDH DRHVIGDPYK TMFLSRLSYD TKESDIEREF TRYGPIERIR VVRNKVTGKS MGYAFVVFER ERDLKVAYKA SAGLMLNGRR IVVDVERGRT VKGWLPRKLG GGLGGRHYTK ERPRRERGSR FRGDSGFRGG YRGGFRKSSG GGSRFGRGPT RSSHSSDYGG RDSSPKRRRY N
Source	E.coli
Target Names	usp101
Protein Names	Recommended name: U1 small nuclear ribonucleoprotein 70 kDa homolog Short name= U1 70K Short name= U1 snRNP 70 kDa homolog Short name= U1-70K Alternative name(s): U1 small nuclear ribonucleoprotein usp101 Short name= U1
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.