



# Recombinant Schizosaccharomyces pombe Cell cycle control protein cwf8 (cwf8)

<b>Product Code</b>	CSB-YP517610SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O14011
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFCSISGETP KEPVISRVSG NVYEKRLIEQ VIRETSKDPV TQQECTLEDL VPVKVPDFVR PRPPSATSLP ALLSLFQEEW DSVALEQFEL RRNLTKQKQ LSTALYSLDA ALRVISRLTK ERDEAREALA KFSDNIGTVS SKTIEVQEVE MGESDDQLKT SLRSTVEKTF QELSSKRKRT KLQPKWATDD AVSQLLQATP STILENLETE STTSFFPSPE NSSFVLCLHK DELLCLDIQS NSTLKIFEQS ALACCWLTSS KIAVATADAI SIFEFPVSSS GLQSVGEPTS SIPIDEKVN LQAHPGSEYL LAASNEKCYI FSLKSQVYNI TVAQHITSLA VHPDGNLFVA GLENLGRFF ETSSGNELTK FGPHSSPVKT LQFGENGYWL VVTTNDDSDI FIWDLRKSEL VQKIPLQTKV AAVSLDITSQ LLVSSDGETL YVHIYVKSSK SWRCMSQTHV SSISNLVWLN ELHQLLFSTS NGAILRLG
<b>Source</b>	Yeast
<b>Target Names</b>	prp19
<b>Protein Names</b>	Recommended name: Cell cycle control protein cwf8
<b>Expression Region</b>	1-488
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.