



Recombinant Schizosaccharomyces pombe Uncharacterized protein UNK4.17 (SPAC2E11.17, SPACUNK4.17)

Product Code	CSB-BP517622SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14082
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MMQNIKDEAY KTYQNVKDEI SSHLPQSAPE PTGPPFKCAI FGCGGINFGT AEGPWNNNSQK LELVLGHLRQ VVALLSPHKS SHERVLKQKA ETEYASSYAN TREFYSADEY LEYLDSPHEE KPDAYIIGIP PETHGSTQKG ADLELAILKR FPDASLFIEK PISAAPVEDA FAVARRLPSE SVISVGYMLR YLKTIVQKAKE IIEKENLKVV STVAKYNSAY IHNSKKFWWI MSESGGPVVE QGTHFCDLR YFGGDVEIDS IKVNRVEWDD PSGKLNAMPV DEKTIPPEER IPRFTAASWK YKGGVGAFT HNITLQGAKY DTCIEVQADG YYLRIVDLYE EPVLYVRSPD SDEEKVYKYP NDDPYYSEFK IFLDAAEGKG DKNLIHSDFL DAVKTYELTK AITEAK
Source	Baculovirus
Target Names	SPAC2E11.17
Protein Names	Recommended name: Uncharacterized protein UNK4.17
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.