



# Recombinant *Schizosaccharomyces pombe* Uncharacterized dipeptidase C3A11.10c (SPAC3A11.10c)

<b>Product Code</b>	CSB-EP517629SXV-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O14124
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	Y LGKGKWIWTP ERKAHFVLTH FPLIDGHNDL PIYLRENYDN RLLNISLEHL PGQTDIFRLR QGHVGGQFWS VFVECPSLDS NSSLSWNRTG EYEAVTQTLQ QIDVVKRMAL YYPKTFSLTD HSGKVKFDL RNHISSMMGI EGLHQIAGSP SILRQFYDLG VRYATLAHNC DNVFADAADV GKRTNKGLSP AGRDIVREM N RLGMIVDLSH TTPETMHQAL DVSVAFAFFS HSSAKGVYDH PRNVPDDVLI RVKETDGVVM VNFYPAFISP HPENATIDTV VEHIMHIANV TGSYRHIGLG GDFDGIDMVP KGLEDVSKYP DLFVKLAERG LSITELADIA GRNVLRVWKT TEDLGHSIHE PPLEWEDDF
<b>Source</b>	E.coli
<b>Target Names</b>	SPAC3A11.10c
<b>Protein Names</b>	Recommended name: Uncharacterized dipeptidase C3A11.10c EC= 3.4.13.19
<b>Expression Region</b>	40-409
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.