



Recombinant Arabidopsis thaliana Transcription factor UPBEAT1 (UPB1)

Product Code	CSB-EP517790DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O22901
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGVTLEGQRK ESIWVLMRRQ RARRALVKKI MIRPRKSVEA SRRPCRAIHR RVKTLKELVP NTKTSEGLDG LFRQTADYIL ALEMKVKVMQ TMVQVLTETN CV
Source	E.coli
Target Names	UPB1
Protein Names	Recommended name: Transcription factor UPBEAT1 Alternative name(s): Basic helix-loop-helix protein 151 Short name= AtbHLH151 Short name= bHLH 151 Transcription factor EN 146 Transcription factor bHLH151 bHLH transcription
Expression Region	1-102
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.