



# Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At4g16470 (PCMP-E12), partial

<b>Product Code</b>	CSB-MP517810DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O23491
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKTLISSLS HKSQCLTICF KLNKLAMVYY TREFQTEASQ TSASGSMFSG NATTILRRML AEKRIGRFQV ENQRKTEKLD KTLKGLCVTG RLKEAVGLLW SSGLQVEPET YAVLLQECKQ RKEYTKGKRI HAQMFVVGFA LNEYLKVKLL ILYALSGDLQ TAGILFRSLK IRDLIPWNAM ISGYVQKGLE QEGLFIYYDM RQNRIVPDQY TFASVFRACS ALDRLEHGKR AHAVMIKRCI KSNIIVDSAL VDMYFKCSSF SDGHRVFDQL STRNVITWTS LISGYGYHGK VSEVLKCFEK MKEEGCRPNP VTFLVLTAC NHGGLVDKGW EHFYSMKRDY GIEPEGQHYA AMVDTLGRAG RLQEAYEFVM KSPCKEHPPV WGSLLGACRI HGNVKLLELA ATKFLELDPT NGGNYVVFAN GYASCGLREA ASKVRKRMEN AGVKKDPGYS QIELQGEVHR FMKDDTSHRL SEKIYKKVHE MTSFFMDIDY YPDGLDSSCP V
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PCMP-E12
<b>Protein Names</b>	Recommended name: Pentatricopeptide repeat-containing protein At4g16470
<b>Expression Region</b>	1-501
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.