



Recombinant *Saccharomyces cerevisiae* Cruciform DNA-recognizing protein 1 (CRP1)

Product Code	CSB-MP518765SVK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	E7NIP0
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain FostersO) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSSELMFNYYT FSWPAGPKDV ILTGTFDDWR GTLPLVKTAK GNFEITMPVK LANKDDTFQF KFIVDGVWCV SDSYKKEHVS EGIENNFLQI TDLVETQEVA GASRIPEAGG LLCGKPPRSA GPPSTSNRKK NKRNNKKRRS KLKKKSTKNN KKSNESSLDDN EEEDGVTGTT TEDVTGTSRE ETPLAEPNTV SKEAPGNFHI LPIDQSADTT QSNGIIGGPG PVLVPNPGEI KEFTEIRDVD ARELNERLNK KEEVPEPVAG PIVESSVTEK SPALPQADDL IVETKEMAHN VQELTPQVEA VTPLINEPEP LPTPEAQISI PESSKVEPVE GSLQPKLVEK RESTEGVSDG SKKVENKAKK DEEVFTLDPI VNKAPKLPLT DEQTAEGRKS PAVSEEKEKK KKQEKGSKEV KRSETSKEKK PSAKEVKKQT VKAPKKQTAS PLSSTEEP KKKTGFFGKL KKLKFK
Source	Mammalian cell
Target Names	CRP1
Protein Names	Recommended name: Cruciform DNA-recognizing protein 1 Cleaved into the following 2 chains: 1. CRP1 short N-terminal subpeptide 2. CRP1 short C-terminal subpeptide
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.