



Recombinant Arabidopsis thaliana B3 domain-containing protein REM-like 2 (At2g24696)

Product Code	CSB-MP518793DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	F4IQL7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLFNVSACDI RDKPMP SNND INKIDQLERV KKVRKNSSQS EAGSSSSGNS SFVALVKASN LKEDALYLPQ DCTSSNGLNK KCRKIFLTDG GDRSWEMDLK FDKNLDSFCI TRGWRHFCDE NGKKLPNQER FVTVRLVPDC LRNKRLYLSR RFLKNNGLGE PKMVTLVGTD GTRILANLLR ESTGRMSLGR GWVDFAKANR LKIGEYFTLE SIWENDSPIL SLYGTNTSKS DKRKRRENFP VACEKEYVST EARNRNEPEK DKNTEEMINQ ASLSENRLVI TLVPEDVKAG MLRLPSHF MK ANGIDKVGKI YMLGINEMEW WWGDLLTRDG IVSVGCGWRY FCESNGVKIG KSFTLECMYK YDTRPVFKFC PKSGK
Source	Mammalian cell
Target Names	At2g24696
Protein Names	Recommended name: B3 domain-containing protein REM-like 2 Alternative name(s): Protein REPRODUCTIVE MERISTEM-like 2
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.